

ABOUT US


Rhein Cosmetic Surgery Solutions, a division of Rhein Medical, Inc., has been developed to address unmet needs in the plastic and reconstructive surgery market. This division specializes in proprietary products designed to accelerate and improve the recovery process for patients undergoing cosmetic surgery. We understand the importance for surgeons to be in more control of post operative surgical results. In addition to achieving optimal results, we understand the importance for patients to be as comfortable as possible throughout the recovery process. The post operative healing process is a critical step in the overall success of a cosmetic surgery procedure. Rhein Cosmetic Surgery Solutions is committed to developing products that equip doctors and patients with tools to help facilitate recovery, reduce pain, swelling and bruising, and thus improve post operative surgical results and overall patient satisfaction.

Our latest initiative, the Rhein Recovery Bra, is a specially designed system that combines compression therapy with cold compress therapy for patients undergoing breast surgery. This system helps to ensure that patients have the proper comfort and care throughout the at-home healing process. Rather than searching for all of the components in multiple stores, this product is available as a comprehensive kit, making pre operative preparation and post operative care available in one convenient and affordable system. For more information and step-by-step instructional videos, please visit www.RheinCosmeticSurgerySolutions.com.

Rhein Cosmetic Surgery Solutions

Post Operative Care For Patients Undergoing Breast Surgery
Accelerate The Recovery Process ~ Recover In Comfort

Introducing The Rhein Recovery Bra & Rhein GEL-O's,
A Convenient System To Provide Compression Therapy &
Hands-Free Cold Therapy For Post Operative Surgical Recovery.



Premium Breathable Stretch Fabric For Superior Comfort.



Easy GEL-O's Insertion For Hands-Free Cold Therapy.



The Rhein Recovery Bra Is Currently Available In:
Small (Up to 32C), Medium (32D-36DDD), Large (Above 38A)



www.RheinCosmeticSurgerySolutions.com

To order, please call
1-800-637-4346



3360 Scherer Drive, Suite B, St. Petersburg FL 33716 • 800-637-4346 •

Tel: 727-209-2244 • Fax: 727-341-8123

Developed In Coordination With Robert D. Rehnke, M.D. & ReCon Surgical.

Patient Information

AHBE



INCLUDES:
Organza Carry Bag, Two Recovery Bras & Two GEL-O's.

**Doctor & Patient
Tested & Approved**



www.RheinCosmeticSurgerySolutions.com

Call 1-800-637-4346
For Order Information



Two Recovery Bras For Compression Therapy



- Light Compression To Help Speed Healing.
- Step-In Design To Avoid Reaching Overhead.
- Front Adjustable Straps To Minimize Reach.
- No Hooks, No Clasps, No Wires.



Two GEL-O's For Cold Compress Therapy



- Reduce Swelling & Bruising.
- Relieve Pain & Discomfort.
- Central Opening Provides No-Cold Zone.
- No Mess, Hands-Free Cold Therapy.

Rhein Recovery Bras have been designed with internal pockets to hold Rhein GEL-O's to aid in applying hands-free cold therapy for post surgical recovery. The pockets prevent direct contact with the skin, making cold therapy comfortable and convenient.

Two bras are included in the system to allow patients to step out of a damp bra and into a dry one after each cold therapy session. Adjustable straps allow patients to loosen to fit GEL-O's into pockets, then tighten for desired compression after cold therapy.

It is recommended to step in and out of the Rhein Recovery Bra to avoid reaching overhead. See photos below.



Step Into The Recovery Bra.



Gently Pull Over Hips.



Adjust For Comfort.

Rhein GEL-O's have been designed to fit into the Rhein Recovery Bra pockets to provide convenient, hands-free cold therapy for post surgical recovery.

Freeze GEL-O's for two hours prior to applying cold therapy. Apply GEL-O's cold therapy in 20 minute intervals, as needed, to reduce swelling, bruising and pain. The central no-cold zone helps prevent cold contact with sensitive area.

It is recommended to insert the GEL-O's into the Rhein Recovery Bra once the bra is in place on the body. Gently adjust GEL-O's and straps for desired comfort. See photos below.



Locate Pocket On The Bra.



Bend At Tabs & Insert Into Pocket.



Adjust For Comfort.

